Hydrofera Blue TRANSFER™

[en] Antibacterial Foam Dressing

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Rx Only

Hydrofera, Hydrofera Blue and Hydrofera Blue Transfer are trademarks of Hydrofera, LLC. ©2022 Hydrofera, LLC. All rights reserved.

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Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

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Instructions for Use [en]

Description

Hydrofera Blue TRANSFER antibacterial foam dressing is a sterile absorptive foam dressing made of polyurethane (PU) foam, methylene blue (less than or equal to 0.00035 g/g), and gentian violet (less than or equal to 0.00035 g/g).

Indications

Hydrofera Blue TRANSFER antibacterial foam dressings are intended as external dressings for use in local management of wounds such as pressure ulcers, donor sites, venous stasis ulcers, arterial ulcers, diabetic ulcers, abrasions, lacerations, superficial burns, post-surgical incisions, and other external wounds inflicted by trauma. Hydrofera Blue TRANSFER antibacterial foam dressings are indicated for moderate to heavily exuding wounds.

Contraindications

Hydrofera Blue TRANSFER antibacterial foam dressings are not indicated for third-degree burns and should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or one of its components.

Application of Hydrofera Blue TRANSFER Antibacterial Foam Dressing

Wound Site Preparation

1) Cleanse the wound in accordance with established procedures.

2) Carefully dry the periwound skin.

Dressing Application

1) Select the appropriate dressing size to ensure the dressing will cover the entire wound (Figure 1).

2) The dressing may be cut in order to optimize placement over wounds on heels, elbows etc.

3) Apply to wound bed directly. Loosely fill shallow-wounds.

4) Cover and secure dressing with a nonocclusive secondary dressing (Figure 2).

5) The dressing can be applied under compression and total contact casting (Figure 3).

Dressing Change and Removal

Hydrofera Blue TRANSFER antibacterial foam dressings can be left in place up to 7 days, depending on the exudate and clinical conditions of the wound. The dressing should be changed if it turns white or upon strikethrough drainage.

If the dressing has retained its color where it is in contact with the wound, the dressing may be left in place for up to seven days.

Precautions and Observations

- This material is not an implant.
- Do not reuse the dressing.
- Infections:** All wounds can be expected to contain some microorganisms; however, if clinical symptoms should develop, appropriate steps as defined by the attending physician should be taken to address that infection. Regular evaluation and cleansing of infected wounds should be a common practice.
- Underlying Causes:** Some wounds, such as leg ulcers and pressure ulcers, can be nonhealing unless steps are taken to correct the underlying pathology. Hydrofera Blue TRANSFER antibacterial foam dressings are designed for local wound management and to be part of the overall management program for these types of wounds.

Note: Hydrofera Blue TRANSFER antibacterial foam dressings have two organic pigments bound to the foam to render it antibacterial. The dressing has been tested in vitro and found to be an antibacterial barrier for the following bacteria and yeast: Methicillin Resistant Staphylococcus aureus (MRSA), Staphylococcus aureus, Staphylococcus epidermidis (coagulase neg.), Pseudomonas aeruginosa, Pseudomonas fluorescens, Escherichia coli, Enterococcus faecalis, Vancomycin Resistant Enterococcus (VRE), Streptococcus pyogenes, Klebsiella pneumonia, Proteus mirabilis, Proteus vulgaris, Enterobacter aerogenes, Bacillus subtilis, Yersinia enterocolitica, Serratia marcescens, Streptococcus pneumoniae, Candida albicans, Candida krusei, and Candida glabrata. All safety and effectiveness testing is on file at Hydrofera LLC.

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Tous les résultats de tests d'innocuité et d'efficacité sont disponibles auprès d'Hydrofera LLC.